

CLIO 8 HW/SW UPGRADE INSTALLATION

➡ THESE INSTRUCTIONS ARE NEEDED TO UPGRADE YOUR CLIO SYSTEM WITH SC-01 SIGNAL CONDITIONER WITH CLIOwin 7 (OR 6.5) SOFTWARE TO THE NEW CLIO 8 USB SYSTEM.

To upgrade your CLIO system follow these instructions:

1 DISCONNECT AND OPEN YOUR SC-01

- 1) Disconnect the SC-01 unit from the PC.
- 2) Open the SC-01 unit unscrewing the 6 upper screws and removing the upper cover.



2 EXCHANGE RS-232 CARD WITH USB CARD

- 1) Inside the SC-01 locate the four screws securing the upper PCB (with RS-232 connector) and unscrew them.
- 2) Gently remove the upper PCB disconnecting it from the lower one.
- 3) Insert the new USB card in the slot and screw it down firmly.



3 CLOSE YOUR UPGRADED SC-01

- 1) Close the SC-01 unit upper cover.
- 2) Now your SC-01 is functionally identical to the new SC-02. **Your CLIO system has a new 10-digit serial number.** Attach the new serial number sticker on the side of your SC-01 in place of the old 8-digit one.
- 3) Reconnect your upgraded SC-01 Signal Conditioner to the PC.



4 READ THE CLIO 8 INSTALLATION INSTRUCTIONS AND START INSTALLING THE CLIO 8 SOFTWARE FROM POINT 3!

Please note:

- 1) You do not need to change your existing PB-4281 PCI card.
- 2) Remember to register your upgraded CLIO system at www.audiomatica.com.